

ANTI-BOVINE CCL2 Polyclonal Antibody



- Description:** Anti-bovine CCL2 polyclonal antibody
- Catalog Number:** PB0092B-100
- Lot Number:** HJ0154EJ
- Quantity:** 100 µg
- Concentration:** 0.8 mg/mL
- Volume:** 125 µL
- Host:** Rabbit
- Alias:** MCP-1
- Purification:** Bovine CCL2 affinity chromatography
- Immunogen:** Recombinant bovine CCL2
- Molecular Weight:** 8.8 kDa (calculated) for bovine CCL2
- Applications:** **Western blot** - 0.1-1 µg/mL
ELISA - 0.5-10 µg/mL
- The recommended concentrations are approximate values.
Optimal dilutions should be determined for each application.
- Formulation:** 1X PBS buffer containing 0.09% NaN₃
- Stability and Storage:** Stable for up to twelve months from date of receipt at 2-8°C.
- Country of Origin:** USA
- Warranty:** Products are warranted by Kingfisher Biotech, Inc. to meet stated product specifications and to conform to label descriptions when used, handled and stored according to instructions. Unless otherwise stated, this warranty is limited to one year from date of sale or expiry date, whichever comes first. Kingfisher Biotech's sole liability for the product is limited to replacement of the product or refund of the purchase price. Kingfisher Biotech products are supplied for research applications. They are not intended for medicinal, diagnostic or therapeutic use. The products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Kingfisher Biotech, Inc.

Quality Control Signature:

A handwritten signature in black ink that reads "Laura J. Kakach". The signature is written in a cursive style and is enclosed in a thin black rectangular border.

REAGENTS FOR VETERINARY RESEARCH

