

# ANTI-SWINE TNF- $\alpha$

## Polyclonal Antibody



- Description:** Anti-swine TNF- $\alpha$  polyclonal antibody
- Catalog Number:** PA0144S-100
- Lot Number:** EK0323EK
- Quantity:** 100  $\mu$ g
- Concentration:** 1.0 mg/mL
- Volume:** 100  $\mu$ L
- Host:** Rabbit
- Purification:** Protein A affinity chromatography
- Immunogen:** Recombinant swine TNF- $\alpha$
- Molecular Weight:** 16.9 kDa (calculated) for swine TNF- $\alpha$
- Applications:** **Western blot** - 0.1-1  $\mu$ g/mL  
The recommended concentrations are approximate values. Optimal dilutions should be determined for each application.
- Cross-reactivity:** Recombinant bovine TNF- $\alpha$  (++)  
(Western blot) Recombinant equine TNF- $\alpha$  (+)
- Formulation:** 1X D-PBS containing 0.01% NaN<sub>3</sub>
- Stability and Storage:** Stable for up to twelve months from date of receipt at 2-8°C.
- Country of Origin:** USA
- Warranty:** Products are warranted by Kingfisher Biotech, Inc. to meet stated product specifications and to conform to label descriptions when used, handled and stored according to instructions. Unless otherwise stated, this warranty is limited to one year from date of sale or expiry date, whichever comes first. Kingfisher Biotech's sole liability for the product is limited to replacement of the product or refund of the purchase price. Kingfisher Biotech products are supplied for research applications. They are not intended for medicinal, diagnostic or therapeutic use. The products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Kingfisher Biotech, Inc.

REAGENTS FOR VETERINARY RESEARCH

