

RECOMBINANT SWINE CCL3L1

KINGFISHER



Catalog Number: RP0075S-025

Quantity: 25 µg

Alias: MIP-1 alpha

Source: Yeast

Formulation/Reconstitution: Lyophilized. Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein.

Stability and Storage: Stable for up to twelve months from date of receipt at -20°C. Stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.

Molecular Weight: 7.8 kDa (calculated)

Purity: >95% as visualized by SDS-PAGE analysis.

Purification: Ion-exchange chromatography

Entrez Gene ID 494459

Amino Acid Sequence: APLGADTPTA CCFSYTSRQL PRKRVADYFE TSSQCCKPGV
IFQTKKGREV CANPEDAWVQ EYISDLELNA (70)

Bioactivity: The biological activity of recombinant swine CCL3L1 was determined in a chemotaxis assay using a Boyden chamber with freshly isolated swine PBMC monocytes using a Percoll gradient enrichment method. Biological activity typically peaks at 100 ng/mL.

Country of Origin: USA

Support: Generation of this material is based upon work supported by the Cooperative State Research, Extension, Education Services (CSREES), USDA under project #0206006 US-Veterinary Immune Reagent Network (www.vetimm.org).

Warranty: Products are warranted by Kingfisher Biotech, Inc. to meet stated product specifications and to conform to label descriptions when used, handled and stored according to instructions. Unless otherwise stated, this warranty is limited to one year from date of sale or expiry date, whichever comes first. Kingfisher Biotech's sole liability for the product is limited to replacement of the product or refund of the purchase price. Kingfisher Biotech products are supplied for research applications. They are not intended for medicinal, diagnostic or therapeutic use. The products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Kingfisher Biotech, Inc.

REAGENTS FOR VETERINARY RESEARCH

