

# RECOMBINANT FELINE SCF

KINGFISHER



**Catalog Number:** RP0150F-005

**Quantity:** 5 µg

**Alias:** Stem Cell Factor; Kit Ligand

**Source:** Yeast

**Formulation/Reconstitution:** Lyophilized. Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein.

**Stability and Storage:** Stable for up to twelve months from date of receipt at -20°C. Stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.

**Molecular Weight:** 27.9 kDa (calculated)

**Purity:** >95% as visualized by SDS-PAGE analysis.

**Purification:** Ion-exchange chromatography

**Entrez Gene ID:** 493937

**Amino Acid Sequence:** GLCRNRVTDD VKDVTKLVAN LPKDYKIALK YVPGMDVLPS  
HCWISVMVEQ LSVSLTDLLD KFSNISEGLS NYSIIDKLVK  
IVDDLVECVE GHSSENVKKS SKSPEPRLFT PEEFFRIFNR  
SIDAFKDLEM VASKTSECVV SSTLSPEKDS RVSVTKPFML  
PPVAASSLRN DSSSSNRKAT NPIEDSSIQW AVMALPACFS  
LVIGFAFGAF YWKKKQPNLT RTVENIQINE EDNEISMLQE  
KEREFAQEV (248)

**Bioactivity:** The biological activity of recombinant feline SCF was measured in a cell proliferation assay using the human TF-1 cell line. The ED50 for this effect is typically 2 - 4 ng/mL.

**Country of Origin:** USA

**Warranty:** Products are warranted by Kingfisher Biotech, Inc. to meet stated product specifications and to conform to label descriptions when used, handled and stored according to instructions. Unless otherwise stated, this warranty is limited to one year from date of sale or expiry date, whichever comes first. Kingfisher Biotech's sole liability for the product is limited to replacement of the product or refund of the purchase price. Kingfisher Biotech products are supplied for research applications. They are not intended for medicinal, diagnostic or therapeutic use. The products may not be resold, modified for resale or used to manufacture commercial products without prior written approval from Kingfisher Biotech, Inc.

REAGENTS FOR VETERINARY RESEARCH

